



first test installations planned for Q2/2024

## How it works

- take a sample
- insert it into the autosampler
- start the analysis
- retrieve the results

## Features

- Table-top device to count bacteria
- Cooled autosampler (10 samples)
- Automated sample processing and analysis
- Simplified user interface with remote control
- Online data provision
- Connection to LIMS interface
- Responsive touchscreen (15.6")
- Dimensions 100 × 40 × 50 cm (l × w × h)

## Testkits\*

For *Legionella pneumophila*

- L.p. SG1
- L.p. SG2-14

contain all consumables for 20 analyses

\*Validated reagents exclude false-positive results due to cross-reactivities with other water pathogens

## Benefits



### Automated

no sample preparation



### Digital

simple data handling



### Fast

time to result ca. 3 h

